<b>ACCUCHEK® Instant</b>
<b>CONTROLS</b>
REF 07869525

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**For in vitro diagnostic use**

## Intended Use

The Accu-Chek Instant control solution is intended for performing control tests on the Accu-Chek Instant and Accu-Chek Instant S blood glucose meters and Accu-Chek Instant test strips. It is intended for self-testing by people with diabetes and for near-patient testing by healthcare professionals.

Testing control solutions with known glucose levels establishes that the operator and the system are performing acceptably. Control results must be within the defined acceptable ranges before valid patient testing is allowed.

## Important Information

**⚠ WARNING**

**Risk of suffocation**

This product contains small parts that can be swallowed. Keep the small parts away from small children and people who might swallow small parts.

**⚠ WARNING**

**Risk of a serious health incident**

Failure to follow testing instructions or test strip storage and handling instructions can lead to an incorrect test result that may lead to improper therapy.

Carefully read and follow the instructions in the User’s Manual and package inserts for the test strips and control solutions.

Do not ingest or inject the control solution or use the control solution as eye drops.

EUH208 – Contains DIAZOLIDINYL UREA. May produce an allergic reaction.\*

\* Regulation (EC) No 1272/2008 (CLP)

**Contents of the pack**

Pack containing 2 × 2,5 mL control solution and package inserts. The Accu-Chek Instant control solution pack contains two control solutions, one for the hypoglycaemic range (control solution 1, gray cap) and one for the hyperglycaemic range (control solution 2, white cap).

**Disponal**

All components of the pack can be discarded in domestic waste. Consult local ordinances as they may vary by country. If you have any questions, contact Roche.

**Control solution storage and handling**

- Refer to the test strip package insert for system operating conditions.
- The printed use by date is valid if the unopened control solution is stored between 2 and 32 °C.
- The control solution does not need to be kept in the refrigerator. Do not freeze the control solution.
- Write the date the control solution bottle was opened on the bottle label. The control solution must be discarded 3 months from the date the control solution bottle was opened (discard date) or on the use by date on the bottle label, whichever comes first.
- Control solution taken directly from the refrigerator must be allowed to adjust to room temperature (without opening the control solution bottle).
- The control solution can stain fabric. Remove stains by washing with soap and water.

### Performing a Control Test

Refer to the meter User’s Manual for instructions on performing a control test.

## Sources of Error

- If the control results are outside the acceptable range, do not use the meter until you solve the problem. Check this list to help solve the problem.
  - Were the test strips or control solutions expired?
  - Did you wipe the tip of the control solution bottle before use?
  - Were the caps on the test strip container and control solution bottle always closed tightly?
  - Was the test strip used immediately after it was removed from the test strip container?
  - Were the test strips and control solutions stored in a cool, dry place?
  - Was the correct control solution level selected when the test was performed?
  - Did you follow the directions?

## Control Testing Intervals

Follow your facility’s policy for control testing intervals.

Control tests should be performed:

- the first time before using the meter for patient blood glucose testing.
- at intervals established by the facility.
- when a new test strip box is opened.
- if the test strip container was left open.
- if test strips were improperly stored.
- if there is a question about a blood glucose result.
- to check the performance of the system.
- if the meter was dropped.

Your facility may require that control tests be successfully performed **after** any of the following occur and **before** patient testing resumes:

- Previous control results were outside the acceptable range.
- Control tests were not performed at the proper interval.

Control results must be within the designated range on the test strip container label or as defined by your facility, before being considered acceptable. Patients can be tested after control tests have been acceptably performed at the proper testing interval.

## Ingredients

Component	Control Solution Level 1 (Low)	Control Solution Level 2 (High)
Glucose	0.06 <span> </span> %	0.30 <span> </span> %
Buffer	4.83 <span> </span> %	4.82 <span> </span> %
Biological salt	3.39 <span> </span> %	3.38 <span> </span> %
Preservative	0.30 <span> </span> %	0.30 <span> </span> %
Non-reactive ingredients	9.99 <span> </span> %	9.97 <span> </span> %
FD&C Blue #1	0.08 <span> </span> %	0.08 <span> </span> %
Percentages in weight/weight		

**Note:** For an explanation of symbols used, refer to the end of this package insert. Visit our website at www.accu-chek.com or contact the local Roche representative for more information.

## Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please refer to the manufacturer and to your national authority.

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