EN Instructions for use

Intended use

The Accu-Chek LinkAssist Plus insertion device is intended for the insertion of the head set of the Accu-Chek FlexLink Plus and Accu-Chek Insight Flex infusion sets into the skin. The insertion device can be reused multiple times. The insertion device is intended to be used by one person only and must not be used by other people.

Before you get started

Read these instructions for use and the instructions for use of the Accu-Chek FlexLink Plus or Accu-Chek Insight Flex infusion set before using the insertion device.

If the insertion device was dropped or exposed to other mechanical stress, check it for damages, e.g. cracks. Do not use the insertion device if it is damaged.

Information regarding any special care you should exercise for the safe and effective use of the device, or to avoid damage to the device that could occur as a result of use, including misuse, is preceded by the following heading:

PRECAUTION

A Release button	C Release lock
B Pretensioning element	D Loading chamber

Inserting the head set into the insertion device

- If the pretensioning element of the insertion device is up, pull the pretensioning element down to the starting position.
- Insert the head set into the insertion device. The blue button of the head set must point towards the blue release lock of the insertion device. The head set is inserted correctly when you hear it **CLICK** into place.
 Tip: The symbol ➞ on the protective film of the adhesive pad (A) indicates the side where you will later connect the transfer set.
- Remove the white protective film (A). Then remove the blue protective film (B).
- Hold the insertion device with the blue side of the pretensioning element facing up. Hold the pretensioning element with your fingers on both sides. Pull the pretensioning element all the way up until you hear it **CLICK** into place. The introducer needle of the head set with the soft cannula snaps out.
- PRECAUTION**
Risk of injury
 After the head set has been inserted into the insertion device, the introducer needle of the head set protrudes.
 Do not point the inserted head set at your face or other people.
- Place the insertion device on the chosen infusion site. Press the insertion device firmly against the skin. The release lock must disappear completely into the insertion device and must no longer be visible.
 Tip: You should be able to comfortably reach the side of the head set where you will later connect the transfer set. This side matches the side of the insertion device without the release button.
- Press the release button to insert the cannula of the head set into the skin.
 If you cannot press the release button, the release lock has not completely disappeared in the insertion device. In this case, press the insertion device more firmly against the skin and try again.
- Remove the insertion device.
 Continue as described in the instructions for use of the infusion set.
- To store the insertion device, pull the pretensioning element down to the starting position.

Storing the insertion device

Store the insertion device with the pretensioning element pulled down to the starting position. Keep the product dry and away from sunlight. Store the product between 5 and 45 °C.

Cleaning the insertion device

Wipe the outside of the insertion device with a standard disinfectant (e.g. 70 % isopropanol) and let it dry. Do not immerse the insertion device in the cleaning agent.

Disposing of the insertion device

The device of the insertion device as plastic waste. For information about correct disposal, contact your local council or authority.

Reporting of serious incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745/EU on medical devices), if, during the use of this device or as a result of its use, a serious adverse incident has occurred, please report it to the manufacturer and to your national authority.

Warranty

This device is covered by the warranty regulations currently valid in your country. Any changes or modifications to the device not expressly approved by Roche could render the warranty null and void.

Customer support

Distributed in the United Kingdom by:

Roche Diabetes Care Limited
 Charles Avenue, Burgess Hill
 West Sussex, RH15 9RY, **United Kingdom**
 Accu-Chek Pump Caroline ¹⁾
 UK Freephone number: 0800 731 22 91
 RO Freephone number: 1 800 88 23 51
¹⁾ calls may be recorded for training purposes
 Some mobile operators may charge for calls to these numbers.
 burgesshill.insulinpumps@roche.com
 www.accu-chek.co.uk
 www.accu-chek.ie

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